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Ion-selective Electrode for Bromide

Introduction

The ISE25Br electrode is a solid-state membrane electrode for direct measurement of bromide ions in aqueous solutions. It is characterised by near theoretical sensitivity, fast response, a high degree of selectivity for bromide ions against a background of other ions, low redox sensitivity, low maintenance level, and sturdy structural properties. The electrode will withstand almost all forms of chemical attack.

Applications

The electrode is primarily intended for measurement of the activity of bromide ions, either by direct measurement or as an indicator electrode for bromide titrations.

In dilute solutions the concentration of free bromide ions equals activity and can be measured directly at a suitable ionic strength, and after proper calibration. The discrepancy between activity and concentration is considerable only in solutions of rather high concentration, where the conversion factor, and the activity coefficient, are not known accurately. Thus, as regards concentrated bromide solutions, for example, 1 M Br^- , it is advisable to dilute the samples, for example 100 times, before measurement. The result is then corrected by multiplying by the dilution factor.

Another problem arises if the sample contains a fairly high background concentration of ions other than Br^- and its accompa-

nying cation. It is recommended that so much of an indifferent salt is added, e.g. KNO_3 , to both standards and samples so that the total ionic strength is almost entirely governed by the KNO_3 addition. When this method does not suffice, and when the sample contains bromide ion complexing agents, calibration curves are useless and other methods must be applied, e.g. the so-called standard addition or standard subtraction methods.

Measuring Hints

For precise measurements, all samples and standards should have the same temperature, preferably ambient temperature. The electrode should not be used above 60°C .

Always rinse the electrode with distilled water between measurements. Careful rinsing is important when the ISE25Br has been exposed to a rather concentrated solution and afterwards has to be used in a very dilute solution. Effective rinsing is achieved by flushing the electrode for about 15 secs. in two lots of distilled water. Dry with a clean, soft tissue to prevent cross contamination.

Measuring Procedure

If the ISE25Br is exposed to strong light, a marked effect on the electrode potential will only be noticeable at concentrations lower than 10^{-3} M Br^- . The light effect is diminished by stirring. Hence, the preferred way of measuring both standard solutions and unknown samples is to pour the solution into a black beaker (part no. 904-515), lower the ISE25Br into the solution, start stirring, and then wait until the meter reading becomes stable. Depending on the concentration and the sample matrix, it will take approximately 60 secs. for the reading to stabilise.

If an ion meter is used, the concentration can be read directly. If a mV reading is obtained, the concentration can be found from the calibration curve.

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Guidelines for Sample Preparation

High ($> 10^{-1}$) and variable ionic strength of sample	
Bromide conc. level	Action
High or mean ($> 10^{-3}$)	Dilute standards and samples 100 times with 0.1 M KNO_3 solution.
Low ($\leq 10^{-3}$)	Reliable direct results not obtainable because ionic strength of standard cannot be properly adjusted. Other methods, e.g. standard addition or subtraction, must be employed.

High ($> 10^{-1}$) and constant ionic strength of sample	
Bromide conc. level	Action
High or mean ($> 10^{-3}$)	Dilute to approx. a 10^{-3} ionic strength level.
Low ($\leq 10^{-3}$)	Add KNO_3 to the standards to obtain same ionic strength as in sample.

Mean (10^{-1} - 10^{-3}) ionic strength of sample	
Bromide conc. level	Action
Mean or low ($\leq 10^{-2}$)	Add KNO_3 to both standards and samples to obtain a high ionic strength level compared with the sample's ionic strength. Ionic strength $> 10^{-2}$.

Low ($<10^{-3}$) ionic strength of sample	
Bromide conc.	Action
Low ($\leq 10^{-3}$)	Standards and samples are measured directly

Ionic strength of sample governed by Br^- and accompanying cation only	
Bromide conc.	Action
High ($> 10^{-2}$)	Dilute sample e.g. 100 times, and multiply direct reading by the same factor.
Mean or low ($\sim 10^{-2}$)	Standards and samples are measured directly.

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
Calibration

If an ion meter is used e.g. PHM250 Ion Analyzer, the electrode can be calibrated to display the bromide ion concentration directly. The PHM250 can perform a calibration on up to 9 calibration points. In other cases, calibration can be based upon mV readings and on a manually drawn calibration curve.

Shown above is an example of a calibration curve obtained with ISE25Br. The measured potentials are drawn on semi-logarithmic paper against the bromide ion concentration of the calibration standards used. To obtain the highest accuracy, the concentration values of the standards used during a calibration should lie on either side of the concentration values of the samples to be measured. Furthermore, the same temperature should be used during calibration and sample measurement.

There is very little need for recalibration owing to the great stability of the ISE25Br. Just how often a calibration should be checked, or a recalibration should be performed depends entirely on the desired accuracy. Due to the fact that the electrode sensitivity (slope) remains unchanged with time, it is normally only necessary to check one point on the calibration curve e.g. by immersing the electrode in a standard solution and comparing the measured potential with a previous measurement.

Special Applications

The ISE25Br is ideal for titration purposes. In fact, the most reliable determinations (± 0.1 % accuracy) of bromide are achieved by potentiometric or coulometric titration, using silver ions as titrant and the ISE25Br as the indicator electrode. Please note that the ISE25Br can be used both in titrations of bromide and in titrations where bromide is the titrant, e.g. in silver determinations.

The ISE25Br has extremely low drift and can, accordingly, be used to advantage for monitoring bromide or for maintaining a constant bromide level in a sample or on-stream, using a technique similar to the well-known pH-stat technique.

Indirect Methods

The **standard addition technique** with one addition is very suitable for analysing samples in which the ionic strength and interferences are moderate to high and variable, or when the sample contains bromide ion complexing agents. With the standard addition technique the electrodes are immersed in the sample and the potential (E) is measured. An aliquot of standard solution containing bromide ions is added to the sample and the potential (E_1) is measured. Although standard addition does not require calibration, the sensitivity of the electrode in the actual measuring range should be known in advance. As a general rule, the addition of standard should double the concentration. For Br^- determination, a potential change between 10 and 30 mV will normally be sufficient.

Upon completion of the measurement, the concentration of Br^- in the sample is calculated as follows:

$$C_{\text{smp}} = \pm \frac{C_{\text{std}} \cdot V_{\text{std}}}{(V_{\text{smp}} + V_{\text{std}}) \cdot 10^{[E_1 - E/S] - V_{\text{smp}}}}$$

Where:

C_{smp} = concentration of Br^- in sample

C_{std} = concentration of Br^- in standard

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V_{smp}	=	sample volume (ml)
V_{std}	=	total volume of standard added (ml)
E	=	potential measured before first addition/ subtraction (mV)
E_1	=	potential measured after addition/ subtraction (mV)
S	=	electrode sensitivity (mV/pC), nominal value = - 59 mV at 25 °C

In addition, the **standard subtraction technique** can be performed. This involves the addition of known amounts of reagents which react with the Br^- in the sample. The same formula (with negative sign) is used for result calculation.

Analate addition techniques can be applied too. During an analate addition measurement, the sample is added to the standard instead of vice versa.

The analate technique is primarily used when the concentration of the ion in the sample is relatively high, thus necessitating a dilution to bring the concentration into the linear area of electrode response. It is also useful in cases where the ionic strength of the sample is very high and the sample volume is very small.

Reference Electrodes

Preferred types:

It is recommended that the ISE25Br is used in combination with the REF251 double junction electrode with 0.1 M KNO_3 in the second salt bridge.

The REF601 Mercury/Mercurous Sulfate Electrode with a salt-bridge solution of saturated K_2SO_4 may also be used.

Salt-bridge solution and liquid-junction:

The above-mentioned salt-bridge solutions are not always recommendable. One requirement of such solutions is that they be able to create a stable liquid-junction potential as near to zero as possible. To achieve this, they must consist of a concentrated salt solution whose anions and cations have approximately the same mobilities. Less concentrated solutions, e.g. 0.1 M KNO_3 can be used, if such a solution is used as supporting electrolyte. Further, bromide, silver and strongly interfering ions must be absent from the salt-bridge, and sparingly soluble precipitates must not be created at the phase boundary. The last problem arises in many organic solvents. In such cases, the REF251 Electrode can be used with a filling of, e.g. LiCH_3COO in alcohol.

Storage

Between measurements and up to one week. Rinse the electrode with distilled water and place it in distilled water.

One week or longer. Rinse with distilled water. Place the electrode in its box and store at room temperature.

Conditioning

When the ISE25Br is taken from dry storage or from storage in distilled water, it is ready for use without any conditioning.

Maintenance

Except for the previously recommended rinsing procedure between measurements, the ISE25Br requires very little maintenance. Only in the event that the electrode has been exposed to strong light for extended periods, or has been lowered into solutions containing strongly interfering ions, will it be necessary to renew the membrane surface by polishing. The polishing procedure is as follows:

1. Polish the membrane by gently pressing it against the emery strip while rotating it for a **few** seconds.
2. The polishing is finished when the membrane obtains its original colour (glaze).

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Limitations

Interferences

The ions that will interfere in bromide determinations with the ISE25Br are listed in Specifications, page 14.

Strongly reducing agents must not be present.

Complexing agents

If bromide ion complexing agents are present, a special measuring technique must be adopted, as previously mentioned, but such agents will not normally influence the electrode's behaviour.

Poisons

The ISE25Br Electrode cannot normally be so badly poisoned by, for example, interfering substances that it cannot be cured by careful polishing of the membrane.

Note: the membrane must never contact wet or dry surfaces of metals less noble than silver.

Light effect

If the ISE25Br is exposed to strong light, a marked effect on the electrode potential will be noticeable only at concentrations lower than 10^{-3} M Br^- . Thus, black beakers (part no. 904-515) must be used for sample measurements.

Stirring diminishes the light effect upon the bromide electrode, although this still exists at low bromide concentrations.

Temperature effect

According to the Nernst equation, the temperature will influence the electrode potential. Except for this unavoidable effect, the ISE25Br is seen to have no hysteresis after its exposure to varying temperatures within the specified range of 0 to 60°C.

Stirring effect

Both the potential of the ISE25Br and the liquid-junction potential of the reference electrode are to a small extent influenced by stirring. To obtain lowest possible light sensitivity and fast hysteresis-free response with almost theoretical sensitivity, it is generally recommended to ensure well-defined stirring conditions during sample measurements. Yet, if stirring is performed by simply shaking the measuring beaker, and if dealing with low bromide content and ionic strength, it is recommended as a precaution to stop stirring before taking the meter reading as soon as it has stabilised. Note that the same methodology must always be used during both standardisation and sample measurements.

Hysteresis

Hysteresis as a result of the "carry over" from one sample to another, because of absorption of ions in the electrode materials or of "dead volumes" in the electrode, is very limited. The previously described rinsing procedure will completely eliminate this problem.


Note: the ISE25Br is in fact constructed so that no "dead volumes" exist.

Detection limit

The detection limit is dependent on the solubility product of AgBr. Accordingly, bromide concentrations lower than 10^{-6} M Br⁻ cannot be determined with the ISE25Br. In fact, if the ISE25Br is lowered into a totally bromide-free solution, it will attain a potential corresponding to a saturated AgBr solution.

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Theory

The ISE25Br is a solid-state membrane electrode. When the membrane is in contact with a solution containing bromide ions, an electrode potential develops across the membrane. This electrode potential is measured against a constant reference potential, using an ion meter. The level of bromide ions, corresponding to the measured potential, E , is described by the Nernst equation:

$$E = E^\circ + RT/zF \cdot \ln a_{\text{Br}^-}$$

where

- E° = standard electrode potential (V)
- R = gas constant ($8.3143 \text{ J} \cdot \text{mol}^{-1} \cdot \text{K}^{-1}$)
- T = temperature (K)
- z = valence (-1)
- F = Faraday constant ($96487 \text{ C} \cdot \text{equiv}^{-1}$)
- a_{Br^-} = activity of the bromide ion

Activity Versus Concentration

As seen from the Nernst equation, it is the bromide ion activity which is measured. The bromide ion activity, a_{Br^-} is related to concentration c_{Br^-} by the following equation :

$$a_{\text{Br}^-} = f_{\text{Br}^-} \cdot c_{\text{Br}^-}$$

Where f_{Br^-} denotes the activity coefficient in very dilute solutions, f_{Br^-} is by definition 1 and concentration equals activity.

For concentrations up to 0.2 M, f_{Br^-} can be determined fairly accurately by means of the Debye-Hückel equation:

$$-\log f_{\text{Br}^-} = A \cdot z^2 \cdot \sqrt{I} / (1 + B \cdot \dot{a} \cdot \sqrt{I})$$

Where I is the total ionic strength of the solution (see the following), and A and B are constants which at 25°C can be set to 0.5115 and 0.3291, respectively. \dot{a} is an ion size parameter which is approximately 3 for Br^- .

By inserting the constants and $z = -1$ (the valency of the bromide ion), the equation can be simplified to:

$$-\log f_{Br^-} = 0.5 \cdot \sqrt{I} / (1 + \sqrt{I})$$

Thus, it can be seen that the ionic strength is an important parameter in measurements with ion-selective electrodes. The ionic strength expresses the concentration of all ions in the solution in terms of both their molarity and their charges. The ionic strength is calculated as follows:

$$I = 1/2 \cdot \sum c_i z_i^2$$

As all ion species present are included in the summation evaluating I , f_{Br^-} will remain almost constant for solutions with a high background concentration of ions. This means that in most cases the divergence between activity and concentration can be overcome by adding an electrolyte which, without interfering, gives a high and constant ionic strength to both the sample and the calibration solutions. In the "Guidelines for Sample Preparation" section, information is provided concerning the addition of electrolyte.

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Accessories

The ISE25Br Electrode part no. E41M001 is delivered complete with protection cap, abrasive strips and operating instructions part no. D31M035.

Reference Electrodes Available

REF251 Reference Electrode with double salt-bridge	945-465
REF601 Reference Electrode with a salt-bridge solution of sat. K_2SO_4	945-470
1 m electrode cable with type 7 plug	A94L116
1 m electrode cable with BNC plug	A94L114
Abrasive strips	991-565

Specifications

Concentration range:	10^{-6} to 10^{-1} M Br
Typical detection limit:	1.5×10^{-6}
pH range:	2 to 12
Typical reproducibility:	± 0.5 mV
Response time:	Typically 60 secs.
Sensitivity (slope):	> 96 %
Temperature range:	0 to 60°C
Checkpoint:	Approx. 0 ± 25 mV at 25°C in 0.01 M solution of KBr in 0.1 M KNO_3 when the REF251 reference electrode is used with 0.1 M KNO_3 in the second salt-bridge.
Sensor type:	Crystalline
Design:	The ISE25Br electrode will resist most reagents, e.g. 25 % HNO_3 , CH_3OH , CH_3CH_2OH and CH_3COCH_3 . However, the following must be avoided: $CHCl_3$, DMF and CH_2Cl_2 .

Dimensions: Diameter: 12 mm
Length: 103 mm

Interference: I⁻, S²⁻, CN⁻

Cations - except those that are precipitated or form complexes with bromide-do not interfere.

Redox effect: Strongly reducing agents must be absent. No redox sensitivity at moderate redox level.

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