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
Reference Electrode

Introduction

The REF321 is a general purpose Ag/AgCl reference electrode, fitted with a screw cap.

Preparation for measurement

1. Remove the protection cap from the electrode and any seals covering the filling hole.
2. Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.
3. Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole. If necessary, refill with KCl•Ag 3 M KCl Solution, saturated with AgCl.

If desired, the concentration of KCl can be changed. However, it is advisable to use a high concentration. Remember always to saturate the solution with AgCl.

For nonaqueous applications such as measurements in acetic acid. Empty the electrode and fill it up again with a saturated KCl solution in acetic acid. Saturate the solution with AgCl.

4. In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the porous plug down.

Maintenance

1. Electrode contamination is a major cause of faulty measurements.
2. The electrode should be rinsed with distilled water after measurements.
3. Check frequently the level of filling solution.
4. In case of deposits which cover the electrode, clean the electrode with:
 - a solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc...
 - KS400 Pepsin in HCl Solution or RENOVO•X Xtra Strong Cleaning Solution: protein deposits (milk, cheese, serums...). Duration of treatment 1 to 2 hrs.

- KS410 Thiourea Solution: for porous plugs contaminated with sulphides or blocked by an AgCl precipitate. Duration of treatment, a few hours until the porous plug turns white.
- RENOVO•N Normal Cleaning Solution : greasy or oily deposits...

The porous plug of the electrode can be cleaned using a fine abrasive paper.

Storage

Between measurements: leave the REF321 in KS110 3M KCl Solution.

Overnight or longer: seal the filling hole with paraffin film or with the electrode clip and fit back in place the protection cap filled with the KS110 3M KCl Solution.

Accessories

| | |
|--|---------|
| KS110 3M KCl Solution, 500 ml | C20C320 |
| KCl•Ag 3M KCl Solution saturated with AgCl, 100 ml | S21M004 |
| ○RENOVO•N Normal Cleaning Solution, 250 ml | S16M001 |
| □RENOVO•X Xtra Strong Cleaning Solution, 250 ml | S16M002 |
| KS400 Pepsin in HCl Solution, 250 ml | C20C370 |
| KS410 Thiourea Solution, 250 ml | C20C380 |
| CL111 Electrode cable with banana plug | A94L111 |

Specifications

| Reference potential when filled with 3M KCl (in mV) | | | | | | |
|---|-----|-----|-----|-----|-----|-----|
| Temperature (°C) | 0° | 10° | 20° | 25° | 30° | 40° |
| vs. std. hydrogen elec. | 225 | 219 | 212 | 208 | 204 | 196 |
| vs. sat. Hg/Hg ₂ Cl ₂ elec. | -35 | -36 | -36 | -36 | -37 | -38 |

Temperature range: 0 to 80 °C

Note: an electrode chain comprising pHG301 or XG100/200 glass electrode and a REF321 filled with 3M KCl will give a zero pH of approx. pH 7.25.

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RADIOMETER ANALYTICAL SAS
D31M070 • Radiometer Analytical SAS • France • 2004-07B