Troubleshooting

First make sure that steps 1 to 4 described in the section "Checking" have been carried out correctly.

| Causes | Solutions |
|---|---|
| <i>Low sensitivity/ slow response</i> | |
| Calibration error | Use fresh, quality buffers. |
| pH glass membrane not hydrated enough or contaminated | Soak the electrode for a few hours in pH 4 buffer. See "Cleaning". |
| Electrode too old | Replace. |
| <i>Drifting potential</i> | |
| No stirring | Ensure reproducible stirring during calibration and measurement. |
| <i>Unstable readings</i> | |
| Liquid junction clogged with crystals | Scrape gently with a needle. |
| Poor connection | Check connection on the meter and that the instrument is properly grounded. |
| Low ionic strength of sample | Shield the beaker and/or add an inert salt. |

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D31M109 • Printed by Radiometer Analytical SAS • France • 2000-06A

RADIOMETER ANALYTICAL SAS
72 rue d'Alsace, 69627 Villeurbanne Cedex, France
E-mail: radiometer@analytical.com Web: www.radiometer-analytical.com
Tel.: +33 (0)4 78 03 38 38 - Fax: +33 (0)4 78 68 88 12

A Hach Company Brand


- Gel -

Combined pH Electrodes Kombinierte pH-Elektroden Electrodes de pH Combinées

Preparing

1.

2.

After long-term storage soak the electrode in 3M KCl (KS110, part no. C20C320) for at least 12 hours to re-hydrate the liquid junction and glass membrane.

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3.

Before performing measurements, check the state of your electrode. See overleaf. For calibration, we recommend the use of Radiometer Analytical pH standards.

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Checking

1.

This electrode is filled with a gel which contains KCl. Therefore there is no need to fill with filling solution.

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2.

Check that the junction opening is not clogged. Scrape gently if necessary.

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4.

Precautions:

- Never regenerate the active tip of the electrode with 1% HF (hydrofluorid acid) solution.
Never use this electrode:
 - In chlorinated media or in organic solvents.
 - For measurements with significant variations in temperature or temperatures over 60°C.
 - If the pH to be measured is lower than 2 or higher than 12.

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Disposal: Ordinary waste, no specific treatment required.

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Storage

Between measurements: lower the electrode in 3M KCl (KS110).

For overnight storage: use 3M KCl (KS110).

For long-term storage (> 72 hours):

1. Rinse the electrode with distilled water.
2. Cover the junction opening with paraffin film. Cover the glass bulb with the protection cap (empty or filled with 3 M KCl).
Never store the electrode in distilled water.
Note: after long-term storage, soak the electrode in 3M KCl (KS110) for approx. 12 hours before use.
3. Place the electrode in its box and store in an upright position (electrode tip down).

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Cleaning

Important: During cleaning, only the glass membrane should be in contact with the liquid. It is also recommended to cover the junction opening with paraffin film.

Clean the electrode weekly with RENOVO•N. Unstable readings can be due to contamination, use one of the appropriate cleaning solutions. Rinse and recalibrate your electrode after each cleaning.

Protein

Use RENOVO•X (S16M002) solution for 3 minutes.

Use KS400 (C20C370) solution of pepsin in HCl for a few hours.

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