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Dear customer,

The acquired test filter set is designed for your own regular checks of your measuring equipment. To make it easier for you, you will find a draft of test records overleaf for your quality surveillance. Transfer the nominal values given in the quality control certificate HDQ178. Make an appropriate number of copies. Should your results exceed the permissible tolerances, please contact your HACH LANGE agent. Any further request should only be addressed to HACH LANGE, Düsseldorf.

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\*) Enter the nominal values and tolerances for your instrument given in the QC certificate HDQ178

Set number: \_\_\_\_\_  
Prüffiltersatz-Nr., Jeu de filtres d'essai no.:

Valid until: \_\_\_\_\_  
Gültig bis: Validité :

Date: \_\_\_\_\_  
Datum: Date de mesure:

Signature: \_\_\_\_\_  
Unterschrift: Signature:

## Test Instruction

for HACH LANGE Instruments DR 6000 / DR 5000 / DR 3900 / DR 3800 / DR 2800 / DR 1900

**Test subject:** Absolute accuracy measured against the golden instrument with test filter set.

**Procedure:** The procedure on your instrument can differ slightly from this description. For details

please see also user manual section Optical Checks > Verification.

Install 10x10mm square cuvette holder. Turn on instrument.

### Straylight test:

1. Select **Single Wavelength** /  $\lambda$ , enter **340**, press **OK**.

2. Press **Zero**.

3. Insert filter **KV450/3** into square cuvette compartment, close the light excluder flap. Press **Read**.

Enter measuring values into the table "actual value", remove filter.

### Photometrical accuracy:

4. Select **Options** /  $\lambda$ , enter **546**, press **OK**.

5. Press **Zero**.

6. Insert filter **NG9/1** into square cuvette compartment, close the light excluder flap. Press **Read**.

Enter measuring values into the table "actual value"/,filter NG9/1", remove filter.

7. Insert filter **NG5/2** into square cuvette compartment, close the light excluder flap. Press **Read**.

Enter measuring values into the table "actual value"/,filter NG5/2", remove filter.

8. Insert filter **NG11/2** into square cuvette compartment, close the light excluder flap. Press **Read**.

Enter measuring values into the table "actual value"/,filter NG11/2", remove filter.

9. Press **Exit** or **Main Menu**.

### Wavelength accuracy of DR 6000 / DR 5000:

10. Select **Wavelength Scan**, **Options** /  $\lambda$ , enter **358** - **362** Step 0.1nm, press **OK**.

11. Press **Zero**.

12. Insert filter **Ho** into square cuvette compartment, close the light excluder flap. Press **Read**. Enter

wavelength of absorbance maxima into the table "actual value"/,filter Ho", remove filter.

13. Select **Options** /  $\lambda$ , enter **806** - **809** Step 0.1nm, press **OK**.

14. Press **Zero**.

15. Insert filter **BG20/2** into square cuvette compartment, close the light excluder flap. Press **Read**.

Enter wavelength of absorbance maxima into the table "actual value"/,filter BG20/2", remove filter.

### Wavelength accuracy of DR 3900 / DR 3800 / DR 2800:

16. Select **System Checks** / **Optical Checks** / **Wavelength Check**. Insert filter **BG20/2** into square

cuvette compartment, close the light excluder flap. Touch **Start**. Enter wavelength of absorbance

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